

FaceTite Mini

For small facial areas

FaceTite Z

For the full face

Compare the Differences



Small areas of the face
(lines between the eyebrows, crow's feet,
fine lines under the eyes, mouth, wrinkles,
nasolabial folds, tear gutters)



Bipolar radio frequency (RF) energy



LED light energy



X



Cordless



X

Forehead, checks, upper lip
jawline and neck



Bipolar radio frequency (RF) energy



LED light energy



Infrared (IR) heat energy



Cordless



Slider Gel

